

APPLICATION

Separation of Six Naturally Occurring Cannabinoids at High pH Using a Kinetex® 2.6 µm EVO C18 Core-Shell Column

Seyed Sadjadi, Laura Snow, and Simon Lomas
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

With the recent legalization of marijuana in various states across the USA, there has been increased focus on the study of medicinal uses of cannabis. There are several main classes of naturally occurring cannabinoids, each with different physical and psychoactive effects. The amount of each cannabinoid present varies depending on the strain of cannabis, thus making accurate quantitation incredibly important. Some strains for example are bred to decrease levels of THC which is psychoactive and increase levels of CBD which has had some success reducing epileptic seizures. Determining the profile for a strain of cannabis can be useful for effective treatment as well as quality control. By utilizing a high performance Kinetex EVO C18 column, excellent resolution and sensitivity are achieved. Moreover, because of the pH 1-12 stability of the Kinetex EVO C18 media, a high pH method can be used; promoting better negative ionization of these acidic compounds.


LC/MS/MS Conditions

Column: Kinetex 2.6 µm EVO C18
Dimensions: 100 x 2.1 mm
Part No.: 00D-4725-AN
Mobile Phase: A: 5 mM Ammonium Bicarbonate (pH 9.5) in Water
 B: Methanol / Acetonitrile (80:20)
Gradient:

Time (min)	% B
0	60
7	95
8	95
8.1	60
10	60

Flow Rate: 0.4 mL/min
Temperature: Ambient
Detection: MS/MS (AB SCIEX 4000 QTRAP®)
Sample:

1. Cannabidiol (CBD)
2. Cannabidiol (CBD)
3. Cannabigerol (CBG)
4. Cannabinol (CBN)
5. Tetrahydrocannabinol (THC)
6. Cannabichromene (CBC)



APPLICATION

Ordering Information Kinetex[®] EVO C18 Core-Shell LC Columns

5 µm Minibore Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298 for 2.1 mm ID

5 µm MidBore [™] Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4633-YO	00D-4633-YO	00F-4633-YO	AJO-9297 for 3.0 mm ID

5 µm Analytical Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJO-9296 for 4.6 mm ID

2.6 µm Minibore Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4725-AN	00B-4725-AN	00D-4725-AN	00F-4725-AN	AJO-9298 for 2.1 mm ID

2.6 µm MidBore [™] Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4725-YO	00D-4725-YO	00F-4725-YO	AJO-9297 for 3.0 mm ID

2.6 µm Analytical Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6		3/pk
00B-4725-E0	00D-4725-E0	00F-4725-E0		AJO-9296 for 4.6 mm ID

1.7 µm Minibore Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298 for 2.1 mm ID

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

Australia
t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China
t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark
t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico
t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands
t: +32 (0)2 503 4015
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico
t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain
t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA** 
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide.
For the distributor in your country, contact Phenomenex USA,
International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Kinetex is a registered trademark and SecurityGuard and MidBore are trademarks of Phenomenex. QTRAP is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX[™] is being used under license.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2015 Phenomenex, Inc. All rights reserved.



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.