


APPLICATIONS

LC/MS/MS Analysis of Ethyl Glucuronide (EtG) and Ethyl Sulfate (EtS) from Urine Using a Luna[®] Omega 1.6 μ m UHPLC Column

Jeff Layne, Brian Rivera, and Simon Lomas
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Brian Rivera
Product Manager

In addition to chromatography, Brian also has a passion for ice cream-making, and enjoys experimenting with bold, new flavors.



Introduction

Ethyl glucuronide (EtG) and ethyl sulfate (EtS) are metabolites of ethanol and are slowly eliminated in the urine. EtG has been measured up to 80 hours after heavy ethanol consumption.

With acceptable limits of quantitation being 500 ng/mL for EtG and 250 ng/mL for EtS¹, an LC column with good polar selectivity is critical since EtG and EtS can elute close to matrix interferences.

In this study, we present a method for fast, accurate and reproducible quantitation and determination of both metabolites using Luna Omega Polar C18, a novel sub-2 μ m UHPLC column with polar selectivity and 100% aqueous stability. The efficiency and selectivity of this novel sub-2 μ m phase allows for an accurate and sensitive method while meeting throughput needs.

Materials and Methods

Reagents and Chemicals

Ethyl- β -D-glucuronide, ethyl sulfate, and deuterated standards (EtS-D₅, EtG-D₅) were purchased from Cerilliant (Cerilliant Corporation, Round Rock, TX, USA).

Sample Preparation

Negative human urine was collected and spiked with ethyl- β -D-glucuronide and ethyl sulfate to prepare a stock standard at 20,000 ng/mL. This standard was serially diluted in urine to make a 5 point calibration curve. Two positive samples were also run.

Each standard and sample was diluted 10-fold with mobile phase A (10 mM ammonium formate in 0.1% formic acid, pH unadjusted) and spiked with internal standard as follows: 100 μ L urine + 885 μ L mobile phase A + 10 μ L EtS-D₅ (10 μ g/mL) + 5 μ L EtG-D₅ (100 μ g/mL)

Experimental Conditions

LC/MS/MS method was performed using a Luna Omega 1.6 μ m Polar C18 column on a Shimadzu[®] Nexera[®] X2 with LC30A solvent module (Shimadzu Scientific Instruments, Columbia, MD, USA) and an upper pressure limit of 1000 bar. MS analysis was performed using a SCIEX 4000 QTRAP[®] LC/MS/MS system (Sciex, Framingham, MA, USA).

LC/MS/MS Conditions

LC Column: Luna Omega 1.6 μ m Polar C18

Dimensions: 100 x 2.1 mm

Part No.: 00D-4748-AN

Mobile Phase: A: 10 mM Ammonium Formate in 0.1% Formic Acid, pH unadjusted (~3.2)
B: 50:50 Methanol/Acetonitrile

| Gradient | Time (min) | % B |
|----------|------------|-----|
| | 0.01 | 0 |
| | 1.00 | 0 |
| | 2.00 | 25 |
| | 2.01 | 95 |
| | 3.00 | 95 |

Flow Rate: 400 μ L/min

Temperature: 40 $^{\circ}$ C

LC System: Shimadzu Nexera X2

Detection: MS/MS

Detection System: SCIEX 4000 QTRAP

Table 1.

Primary and Secondary EtG/EtS MRM Transitions

| ID | Q1 | Q3 |
|-------------------|-------|------|
| EtG 1 (Primary) | 221.2 | 75.0 |
| EtS 1 (Primary) | 124.9 | 80.1 |
| EtG_D5 | 226.1 | 85.1 |
| EtS_D5 | 130.0 | 80.1 |
| EtG 2 (Secondary) | 221.2 | 85.1 |
| EtS 2 (Secondary) | 124.9 | 97.0 |



Result and Discussion

The use of the Luna® Omega 1.6 µm Polar C18 column allowed for the fast elution of EtG and EtS in less than 2 min (**Fig. 1**), with total run times including column cleaning less than 5 minutes. This fast separation allows for multiplexing techniques to increase throughput. Note that the retention of EtS is ~0.9 minutes, and EtG is ~3 minutes. This ensures that low responding EtG elutes away from major a matrix component that responds in the EtS channel, which is a known cause for signal suppression.

Calibration curves were generated over a concentration of 50 ng/mL to 5000 ng/mL. A quadratic regression was used to determine relative response versus concentration using peak area of EtG and EtS/peak area of IS, with 1/X weighting factor), Correlation coefficient for EtS calibration curves are 0.997 and 0.9994 for primary and secondary MRM's, respectively. Correlation coefficient for EtG calibration curves are 0.9998 and 0.9998 for primary and secondary MRM's, respectively.

Quantitative data for EtS and EtG standards are summarized in **Table 2** and **Table 3**. Sample concentrations were calculated using the primary MRM channels, and results are summarized in **Table 4**.

Table 2.
EtS Quantitation

| Sample ID | Analyte Concentration (ng/mL) | Calculated Concentration (ng/mL) | Accuracy (%) |
|----------------------|-------------------------------|----------------------------------|--------------|
| Primary Standard 1 | 50 | 46 | 91.2 |
| Primary Standard 2 | 100 | 109 | 109.0 |
| Primary Standard 3 | 500 | 506 | 101.0 |
| Primary Standard 4 | 1000 | 987 | 98.7 |
| Primary Standard 5 | 5000 | 5000 | 100.0 |
| Secondary Standard 1 | 50 | 54 | 108.0 |
| Secondary Standard 2 | 100 | 91 | 91.4 |
| Secondary Standard 3 | 500 | 502 | 100.0 |
| Secondary Standard 4 | 1000 | 1000 | 100.0 |
| Secondary Standard 5 | 5000 | 5000 | 100.0 |

Table 3.
EtG Quantitation

| Sample | Analyte Concentration (ng/mL) | Calculated Concentration (ng/mL) | Accuracy (%) |
|----------------------|-------------------------------|----------------------------------|--------------|
| Primary Standard 1 | 50 | 54 | 107.0 |
| Primary Standard 2 | 100 | 94 | 94.0 |
| Primary Standard 3 | 500 | 482 | 96.4 |
| Primary Standard 4 | 1000 | 1030 | 103.0 |
| Primary Standard 5 | 5000 | 4990 | 99.8 |
| Secondary Standard 1 | 50 | 54 | 108.0 |
| Secondary Standard 2 | 100 | 91 | 91.4 |
| Secondary Standard 3 | 500 | 502 | 100.0 |
| Secondary Standard 4 | 1000 | 1000 | 100.0 |
| Secondary Standard 5 | 5000 | 5000 | 100.0 |

Table 4.
Sample Results

| Sample Name | Analyte Peak Name | Calculated Concentration (ng/mL) | Analyte Peak Area (counts) | IS Peak Area (counts) |
|-------------|-------------------|----------------------------------|----------------------------|-----------------------|
| Sample 1 | ETS-1 | 763 | 6.75E+04 | 6.76E+04 |
| Sample 2 | ETS-1 | 405 | 4.14E+04 | 7.50E+04 |
| Sample 1 | ETG-1 | 5000 | 4.69E+04 | 6.35E+04 |
| Sample 2 | ETG-1 | 2230 | 2.66E+04 | 5.93E+04 |

Figure 1.
Extracted Ion Chromatograms for EtG/EtS and their Deuterated Internal Standards (EtG-D₃ and EtS-D₃), 500ng/mL




Figure 2.
Extracted Ion Chromatograms for Positive Sample 1





Figure 3.
Extracted Ion Chromatograms for EtG and EtS, 100 ng/mL



Conclusions

A method using a Luna Omega 1.6µm Polar C18, a new UHPLC column with excellent polar selectivity, can be used for accurate and quantitative analysis for ethanol metabolites ethyl glucuronide and ethyl sulfate. The method shows good linearity and accuracy from the concentration ranges of 50 to 5000 ng/mL.

Another method was also presented here as proof of concept for mobile phase optimization. This method also showed good separation for both EtG and EtS, with run times of less than 5 minutes, and could be used for further studies for increasing sensitivity.

References:

1. Jatlow, Peter I. et al. "Ethylglucuronide and Ethyl Sulfate Assays in Clinical Trials, Interpretation and Limitations: Results of a Dose Ranging Alcohol Challenge Study and Two Clinical Trials." Alcoholism, clinical and experimental research 38.7 (2014): 2056–2065. PMC. Web. 4 May 2016.

Further Method Development

Another method was implemented using different mobile phase conditions, to compensate for differences in instrumentation and method requirements.

LC/MS/MS method was performed on an Agilent® 1260 (Agilent Technologies, Santa Clara, CA, USA) with an upper pressure limit of 600 bar. MS analysis was performed using a SCIEX API 4000™ MS/MS.

LC/MS/MS Conditions

- LC Column:** Luna Omega 1.6 µm Polar C18
- Dimensions:** 100 x 2.1 mm
- Part No.:** 00D-4742-AN
- Mobile Phase:** A: 10 mM Ammonium Formate pH 3
Formic Acid, pH unadjusted (~3.2)
B: Acetonitrile with 0.1% Formic Acid

| Gradient: | Time (min) | % B |
|-----------|------------|-----|
| | 0 | 0 |
| | 1 | 50 |
| | 1.1 | 0 |
| | 5 | 0 |

- Flow Rate:** 300 µL/min
- Temperature:** 30°C
- LC System:** Agilent 1260
- Detection System:** SCIEX API 4000



APPLICATIONS

Ordering Information

Luna[®] Omega

| Phases | 1.6 μ m Minibore Columns (mm) | | | | SecurityGuard [™] ULTRA Cartridges [†] |
|-----------|-----------------------------------|-------------|-------------|-------------|---|
| | 30 x 2.1 | 50 x 2.1 | 100 x 2.1 | 150 x 2.1 | 3/pk |
| Polar C18 | 00A-4748-AN | 00B-4748-AN | 00D-4748-AN | 00F-4748-AN | AJO-9505 |
| C18 | 00A-4742-AN | 00B-4742-AN | 00D-4742-AN | 00F-4742-AN | AJO-9502 |

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000



If Luna analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Luna is a registered trademark and SecurityGuard is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. QTRAP is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX[™] is being used under license. Shimadzu and Nexera are registered trademarks of Shimadzu Corporation.

Disclaimer

Phenomenex is not affiliated with Agilent Technologies. Comparative separations may not be representative of all applications

© 2016 Phenomenex, Inc. All rights reserved.