

## Improved Clean Up and Recovery of Pharmaceutical Compounds From Plasma using Strata<sup>™</sup>-X Solid Phase Extraction (SPE) vs. Traditional Liquid-Liquid Extraction Methods


Deborah Jarrett  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Liquid-liquid extraction (LLE) is perhaps the most established clean up technique used in the chromatography field. Although it has been used for years, newer techniques with improved specificity towards particular analytes have allowed analysts to improve recovery and reproducibility of their samples. This work explores the benefits of solid phase extraction (SPE) as compared to LLE in a pharmaceutical setting. It was found that SPE provides cleaner extracts, higher recoveries, and better reproducibility which can greatly improve results.


### Introduction

Diclofenac (**Figure 1**) is a slightly acidic ( $pK_a = 4.0$ ) non-steroidal anti-inflammatory drug (NSAID) that has been used as a post operative pain reliever in adult and pediatric patients. As a pain reliever, diclofenac is purported to act via the cyclooxygenase (COX) pathway by inhibition of prostaglandins. However, it also acts as a NSAID that inhibits the lipooxygenase pathway via a process that diminishes the formation of pro-inflammatory hormones. The subsequent quantization of diclofenac from the small volumes of biological matrices, such as plasma, has been of significant concern. Therefore, this article explores two popular extraction methods, solid phase extraction (SPE) and liquid-liquid extraction (LLE), for the isolation of diclofenac from plasma, using a water matrix as the control.

**Figure 1.**  
Structure of Internal Standard Flurbiprofen and Diclofenac



Flurbiprofen ( $pK_a = 4.2$ )



Diclofenac ( $pK_a = 4.0$ )

### Materials and Methods

The plasma pre-treatment step was the same for SPE and LLE and was comprised of filtration through a gauze cloth. Afterwards, 500  $\mu$ L of diclofenac, which was dissolved in 5 % Methanol, was added to 500  $\mu$ L of plasma, and the solution mixture was then acidified with 600  $\mu$ L of 1M Phosphoric acid.

### Solid Phase Extraction

The pre-treated plasma samples were further cleaned up and concentrated using SPE.

**Cartridge:** Strata-X 30 mg/ 1 mL  
**Part No.:** 8B-S100-TAK  
**Condition:** 1 mL Methanol  
**Equilibrate:** 2 mL Water  
**Load:** 1.6 mL Pre-treated plasma  
**Wash:** 1 mL 5 % Methanol  
**Dry:** 1 minute under vacuum at 10 inches Hg  
**Elute:** 1 mL Methanol  
**Dry down:** Dry down @ 53 °C under a stream of nitrogen for 20 minutes  
**Reconstitute:** Reconstitute in 500  $\mu$ L of mobile phase

### Liquid-Liquid Extraction


After pre-treatment, 5 mL of Hexane:IPA (95:5) was added to the pre-treated solution, which was followed by 1 minute of vortexing, and 10 minutes of centrifugation at 2,000 rpm. Subsequently, 4 mL of the top organic layer was transferred to a clean glass centrifuge tube and then evaporated to dryness under a stream of nitrogen at 53 °C for 20 minutes.

# TN-0044

## APPLICATIONS

### LC/UV

Chromatogram of Diclofenac and IS after SPE extraction from a plasma matrix.



Diclofenac spiked plasma sample (50 µg/mL) after extraction with Strata™-X. Flurbiprofen (IS) was added post-extraction at a concentration of 160 µg/mL. Note: the flurbiprofen was added post blow down, which is also post-extraction

**Column:** Kinetex® 2.6 µm C18  
**Dimensions:** 75 x 4.6 mm  
**Part No.:** 00C-4462-E0  
**Mobile Phase:** Water (adjusted to pH 3.3 with phosphoric acid): Methanol (63:37)  
**Gradient:** Isocratic hold for 5 min  
**Flow Rate:** 1.5 mL/min  
**Temperature:** 30 °C  
**Detection:** UV @ 254 nm  
**Sample:** 1. Flurbiprofen (IS)  
 2. Diclofenac

### Results and Discussion

Traditionally, liquid-liquid extraction (LLE) has been used as the standard clean up procedure for a variety of pharmaceutical samples. Although the technique has been used for many decades it could stand to be improved upon. While performing this work, it was discovered that solid phase extraction provided many be-

nefits over this traditional method including improved recoveries, time and solvent savings, as well as more consistent results.

Plasma samples spiked with diclofenac were used to compare extraction efficiencies of both SPE (using a polymer-based SPE sorbent, Strata-X) and LLE. It was found that SPE on the Strata-X sorbent yields approximately 86 % absolute recovery of 15 µg/mL of diclofenac in the plasma matrix as opposed to 46 % for LLE (Table 1). This accounts for almost a two-fold decrease in recovery when using LLE. Furthermore, while this procedure for LLE involved one extraction step, as opposed to multiple extraction steps with the hexane:IPA mixture, in order to have obtained a greater yield using LLE, a greater amount of solvent would have been required for LLE. This would not have only increased the time to obtain a higher extraction yield, but would also increase the total time required for the evaporation of the solvent.

**Table 1.**  
% Absolute Recovery for Diclofenac


|     | Spiked concentration | Diclofenac | Mean % RSD |
|-----|----------------------|------------|------------|
| SPE | 15 µg/mL             | 86 % (n=4) | 10         |
| LLE | 15 µg/mL             | 46 % (n=4) | 35         |

In addition to SPE providing a greater absolute percent recovery by two-fold over LLE, the Strata-X sorbent procedure shows less variability between cartridges, with a mean % RSD of 10 % while the LLE procedure produced more variability with a mean % RSD of 35 times as much, at 35 % (Figures 2 and 3). According to the % RSD values for SPE and LLE, SPE is more precise and reproducible than LLE for the extraction of pharmaceutical compounds.


# TN-0044

## APPLICATIONS

**Figure 2.**  
Diclofenac Extracted Reference Curve: Solid Phase Extraction in Plasma Matrix



**Figure 3.**  
Diclofenac Extracted Reference Curve: Liquid-Liquid Extraction in Plasma Matrix




### Conclusion


The Strata™-X SPE sorbent utilizes the chemical bonding properties of a pyrrolidone ligand, which makes it a prime target for the retention of hydrophobic neutral compounds, while also retaining basic and acidic compounds under strong organic wash conditions. With the presence of the phenyl ring in the pyrrolidone ligand, pi-pi interactions are prevalent for compounds with aromatic structures, and acidic and basic compounds are able to bind to the nitrogen of the secondary amine and the carbonyl carbon, respectively (**Figure 4**).

**Figure 4.**  
Binding Properties of Strata-X


Pi-Pi Bonding



Hydrogen Bonding,  
Dipole-Dipole Interactions



Hydrophobic Interactions



◆ Analyte-sorbent interaction site highlighted in blue.

LLE instead utilizes two immiscible solvents that compete for interaction with the analyte of interest. While LLE has a universal choice of extraction, SPE poses many advantages over LLE for the extraction of diclofenac, via the existence of multiple binding sites on the Strata-X sorbent. Consequently, this data shows that SPE provides greater absolute recovery of diclofenac when compared to LLE; and is less time-intensive, consumes less solvent than traditional LLE procedures, and provides better reproducibility, thereby demonstrating that the extraction method of choice for pharmaceuticals, such as diclofenac, is SPE.

# TN-0044

# APPLICATIONS

- Australia**  
t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com
- Austria**  
t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com
- Belgium**  
t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
beinfo@phenomenex.com
- Canada**  
t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com
- Denmark**  
t: 4824 8048  
f: 4810 6265  
nordicinfo@phenomenex.com
- Finland**  
t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com
- France**  
t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com
- Germany**  
t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com
- India**  
t: 040-3012 2400  
f: 040-3012 2411  
indiainfo@phenomenex.com
- Ireland**  
t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com
- Italy**  
t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com
- Luxembourg**  
t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com
- Mexico**  
t: 001-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com
- The Netherlands**  
t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com
- New Zealand**  
t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com
- Norway**  
t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com
- Puerto Rico**  
t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com
- United Kingdom**  
t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com
- All other countries:**  
**Corporate Office USA**   
t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

**www.phenomenex.com**  
Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

## Ordering Information

### Strata™-X

| Sorbent Mass | Part No.     | Unit           |
|--------------|--------------|----------------|
| 10 mg        | 8B-S100-AAK  | 1 mL (100/box) |
| 30 mg        | 8B-S100-TAK* | 1 mL (100/box) |
| 30 mg        | 8B-S100-TBJ  | 3 mL (50/box)  |
| 60 mg        | 8B-S100-UBJ* | 3 mL (50/box)  |
| 100 mg       | 8B-S100-EBJ  | 3 mL (50/box)  |
| 100 mg       | 8B-S100-ECH  | 6 mL (30/box)  |
| 200 mg       | 8B-S100-FBJ† | 3 mL (50/box)  |
| 200 mg       | 8B-S100-FCH  | 6 mL (30/box)  |
| 500 mg       | 8B-S100-HBJ  | 3 mL (50/box)  |
| 500 mg       | 8B-S100-HCH  | 6 mL (30/box)  |

### Giga™ Tube

|        |              |                |
|--------|--------------|----------------|
| 200 mg | 8B-S100-FDG  | 12 mL (20/box) |
| 500 mg | 8B-S100-HDG  | 12 mL (20/box) |
| 1 g    | 8B-S100-JDG† | 12 mL (20/box) |
| 1 g    | 8B-S100-JEG  | 20 mL (20/box) |
| 2 g    | 8B-S100-KEG  | 20 mL (20/box) |
| 5 g    | 8B-S100-LFF  | 60 mL (16/box) |

### 96-Well Plate

|       |             |              |
|-------|-------------|--------------|
| 10 mg | 8E-S100-AGB | 2 Plates/Box |
| 30 mg | 8E-S100-TGB | 2 Plates/Box |
| 60 mg | 8E-S100-UGB | 2 Plates/Box |

\* Tab-less tube

† Available in Teflon® tubes

Additional sizes available. Contact your Phenomenex Sample Preparation Specialist for additional information.

## Kinetex® Core-Shell HPLC/UHPLC Columns

### 1.7 µm Minibore Columns (mm)

|            | 50 x 2.1    | 100 x 2.1   | 150 x 2.1   |
|------------|-------------|-------------|-------------|
| <b>C18</b> | 00B-4475-AN | 00D-4475-AN | 00F-4475-AN |

### 2.6 µm Minibore Columns (mm)

|            | 50 x 2.1    | 75 x 2.1    | 100 x 2.1   | 150 x 2.1   |
|------------|-------------|-------------|-------------|-------------|
| <b>C18</b> | 00B-4462-AN | 00C-4462-AN | 00D-4462-AN | 00F-4462-AN |

### 2.6 µm Solvent Saver MidBore™ Columns (mm)

|            | 50 x 3.0    | 75 x 3.0    | 100 x 3.0   | 150 x 3.0   |
|------------|-------------|-------------|-------------|-------------|
| <b>C18</b> | 00B-4462-Y0 | 00C-4462-Y0 | 00D-4462-Y0 | 00F-4462-Y0 |

### 2.6 µm Analytical Columns (mm)

|            | 50 x 4.6    | 75 x 4.6    | 100 x 4.6   | 150 x 4.6   |
|------------|-------------|-------------|-------------|-------------|
| <b>C18</b> | 00B-4462-E0 | 00C-4462-E0 | 00D-4462-E0 | 00F-4462-E0 |

More dimensions and phases available, please inquire.



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to another product of the same phase and comparable dimensions, return the product with comparative data within 45 days for a FULL REFUND.

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

### Trademarks

Kinetex is a registered trademark of Phenomenex, Inc. in the United States, European Union, and other jurisdictions. Strata-X, Giga, and MidBore are trademarks of Phenomenex, Inc. Teflon is a registered trademark of E.I. du Pont de Nemours and Co.

### Disclaimer

Comparative separations may not be representative of all applications

Strata-X is patented by Phenomenex, Inc. U.S. Patent No. 7,119,145

© 2011 Phenomenex, Inc. All rights reserved.