


# APPLICATIONS

## Optimization of Mobile Phase Composition for the Analysis of Synthetic Oligonucleotides using a Clarity<sup>®</sup> Oligo-XT Column

Brian Rivera and Matthew Trass  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



**Brian Rivera**  
Product Manager

*In addition to chromatography, Brian also has a passion for ice cream-making, and enjoys experimenting with bold, new flavors.*



Mobile phase composition for the analysis of three different synthetic oligonucleotides models (Bridged Nucleic Acid, Antisense RNA, and DNA phosphorothioate) was studied and optimized using a Clarity 2.6 $\mu$ m Oligo-XT C18 column.

### Introduction

The characterization of OGNs is important in the drug development process and one common method involves liquid chromatography with detection by either UV (260 nm) or mass spectrometry (LC/MS or LC/MS/MS).

Mobile phase composition is important for the electrospray desorption efficiency and traditional mobile phase composition uses a combination of the alkylamine ion-pairing agent triethylamine (TEA) and the acidic modifier hexafluoroisopropanol (HFIP). Although optimal for LC/MS, TEA and HFIP typically also give better selectivity for OGNs even when UV-Vis is used as the detection method.

Previous studies have shown that a balance of TEA and HFIP are required to ensure both good peak shape and good electrospray ionization efficiency (1). Additionally, work by Chen et al (2) have shown improvements in electrospray desorption efficiency using other alkylamine ion-pairing reagents, including diisopropylethylamine (DIEA).

In this application note, we investigate the use of two mobile phase compositions - TEA HFIP and DIEA HFIP - with Clarity Oligo-XT, a novel core-shell material optimized for synthetic oligonucleotide applications. Three OGN models were investigated: Bridged Nucleic Acid, Antisense RNA, and DNA phosphorothioate. These OGN models were used to be representative of modifications commonly used for therapeutic applications. Methods were developed for optimal peak shape and separation of n-1 failure sequences.

### Materials and Methods

#### Reagents and Chemicals

HFIP, TEA, and DIEA were purchased from Sigma Aldrich (St. Louis, MO, USA).

Crude, desalted OGNs (BNA, 2-MOE Gapmer and Phosphorothioate) were purchased from Integrated DNA Technologies (Coralville, IA, USA). Sequences, along with modifications, are indicated in **Table 1**, above. Samples were prepared at 30 $\mu$ g/mL in nuclease-free water.

**Table 1.**  
OGNs Sequences and Modifications

Oligo Type	Length	Sequence
Bridged Nucleic Acid	19mer	5'-TA/BNA-A/TA/BNA-meC/GT/BNA-T/TA/BNA-T/AC/BNA-G/CC/BNA-C/A-3'
2'-MOE Gapmer	20mer	5'-mC mG mA mC mU A T A C G C G C A A mU mA msU mG mG -3'
DNA Phosphorothioate	19mer	5'-ACT*G*A*C*T*G*A*C*T*G*A*C*G*T*A*C*T-3

### Experimental Conditions

LC method was performed using a Clarity 2.6 $\mu$ m Oligo-XT C18 column on an Agilent<sup>®</sup> 1200 (Agilent Technologies, Santa Clara, CA USA) with an upper pressure limit of 600 bar. Because the intent was to assess chromatographic performance with various mobile phases, UV-Vis was used as the detection method.

HPLC conditions were developed for the two methods using DIEA or TEA as ion-pair. Various gradient profiles and mobile phase compositions were investigated with only the optimized methods shown below.

### HPLC Conditions

#### HPLC Method 1: DIEA and HFIP

**LC Column:** Clarity 2.6 $\mu$ m Oligo-XT  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4746-AN  
**Mobile Phase:** A: 50 mM HFIP & 5 mM DIEA in Water  
B: 50 mM HFIP & 5 mM DIEA in Acetonitrile  
**Gradient:** 5 – 20 % B in 15 min to 95 % wash for 5 min  
**Flow Rate:** 0.3 mL/min  
**Temperature:** 60 °C  
**LC System:** Agilent 1200  
**Detection:** UV-Vis @ 260 nm


#### HPLC Method 2: TEA and HFIP

**LC Column:** Clarity 2.6 $\mu$ m Oligo-XT  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4746-AN  
**Mobile Phase:** A: 100 mM HFIP & 4 mM TEA in Water  
B: 100 mM HFIP & 4 mM TEA in Methanol  
**Gradient:** 10 – 25 % B in 15 min to 95 % wash for 5 min  
**Flow Rate:** 0.3 mL/min  
**Temperature:** 60 °C  
**LC System:** Agilent 1200  
**Detection:** UV-Vis @ 260 nm




**Results**


**Figure 1.**  
Chromatogram for BNA, TEA-HFIP Mobile Phase




**Figure 2.**  
Chromatogram for BNA, DIEA-HFIP Mobile Phase




**Figure 3.**  
Chromatogram for DNA Phosphorothioate, TEA-HFIP Mobile Phase




**Figure 4.**  
Chromatogram for DNA Phosphorothioate, DIEA-HFIP Mobile Phase



**Figure 5.**  
Chromatogram for 2'-MOE Gapmer, TEA-HFIP Mobile Phase



**Figure 6.**  
Chromatogram for 2'-MOE Gapmer, DIEA-HFIP Mobile Phase



## Discussion

The scope of method development was optimization of mobile phase composition and gradient profiles for good peak shape and separation on n-1 failure sequences. With all three modified OGNs, the DIEA ion-pairing was optimal.

As seen in **Figures 1 & 2**, TEA and DIEA were relatively similar for BNA. Good peak shape and n-1 separation were observed. However, for the 19mer DNA Phosphorothioate, only DIEA gave partial separation of n-1/n+1 failure sequences (**Figure 4**). Although DIEA did not give n-1 resolution between 2'-MOE Gapmer, it was able to give better peak shape and slight shouldering of the failure sequence (**Figure 6**).

It is important to note that DIEA, being the more hydrophobic alkylamine ion pair, required acetonitrile as the strong solvent since methanol was unable to elute the strongly retaining OGNs. Conversely, TEA performed better with the methanol as the strong solvent. Care should be taken to use the proper ion-pair and strong solvent for good retention and capacity factor.

Finally, only chromatography was assessed; further experiments would need to be determined to investigate the effect of DIEA and TEA in regards to electrospray desorption efficiency.

## Conclusions

Mobile phase composition is important when developing methods for characterization of OGNs. Two ion-pairing agents, TEA and DIEA, were investigated, with DIEA giving overall better peak shape and selectivity for failure sequences. Further method development considerations include investigation on the effect of HFIP and DIEA/TEA for electrospray desorption.

## References:

1. Rudge, James et al. "Preparation and LC/MS Analysis of Oligonucleotide Therapeutics from Biological Matrices." *Chromatography Today*. 1 Mar. 2011: 16-20. Print.
2. Chen, Buyun, and Michael G. Bartlett. "Evaluation of Mobile Phase Composition for Enhancing Sensitivity of Targeted Quantification of Oligonucleotides Using Ultra-high Performance Liquid Chromatography and Mass Spectrometry: Application to Phosphorothioate Deoxyribonucleic Acid." *Journal of Chromatography A* 1288 (2013): 73-81

## Ordering Information

### Clarity®

1.7 µm Minibore Columns (mm)			SecurityGuard™ ULTRA Cartridges†
Phase	50 x 2.1	100 x 2.1	3/pk
Oligo-XT	00B-4747-AN	00D-4747-AN	AJ0-9515
For 2.1 mm ID			

2.6 µm Minibore and Analytical Columns (mm)					SecurityGuard ULTRA Cartridges†	
Phase	50 x 2.1	100 x 2.1	50 x 4.6	100 x 4.6	3/pk	3/pk
Oligo-XT	00B-4746-AN	00D-4746-AN	00B-4746-E0	00D-4746-E0	AJ0-9515	AJ0-9514
					For 2.1 mm ID	For 4.6 mm ID

5 µm Minibore and Analytical Columns (mm)			SecurityGuard ULTRA Cartridges†	
Phase	50 x 2.1	150 x 4.6	3/pk	3/pk
Oligo-XT	00B-4745-AN	00F-4745-E0	AJ0-9515	AJ0-9514
			For 2.1 mm ID	For 4.6 mm ID

5 µm Semi-Preparative Columns (mm)				SecurityGuard SemiPrep Cartridges*
Phase	50 x 10	100 x 10	150 x 10	3/pk
Oligo-XT	00B-4745-N0	00D-4745-N0	00F-4745-N0	AJ0-9516
For 10 mm ID				

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges**	SecurityGuard PREP Cartridges***
Phase	100 x 21.2	150 x 21.2	250 x 21.2	150 x 30	/ea	/ea
Oligo-XT	00D-4745-P0-AX	00F-4745-P0-AX	00G-4745-P0-AX	00F-4745-U0-AX	AJ0-9517	AJ0-9518
					For 21.2 mm ID	For 30 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000

\* SemiPrep SecurityGuard Cartridges require holder, Part No.: AJ0-9281

\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8223

\*\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8277



If Clarity Oligo-XT analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.



# APPLICATIONS

## Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

## China

t: +86 (0)20 2282-6668  
f: +86 (0)20 2809-8130  
chinainfo@phenomenex.com

## Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

## India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

## Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

## Trademarks

Phenomenex and Clarity are registered trademarks of Phenomenex. Axia and SecurityGuard are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

## Disclaimer

Phenomenex is not affiliated with Agilent Technologies.

Comparative separations may not be representative of all applications

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Clarity Oligo-XT is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2016 Phenomenex, Inc. All rights reserved.